

DEZISEPT Gel

Gel Disinfectant for Hygienic Disinfection of Hands with Panthenol

COMPOSITION

Active substances:

- ethanol (CAS: 64-17-5) 75.67%;
- hydrogen peroxide (CAS: 7722-84-1) 0.15%

CHARACTERISTICS

DEZISEPT Gel is a disinfectant of green colour, pleasant odour, which is applied:

* for hygienic disinfection of hands of medical and pharmaceutical staff; in the food industry and catering; and is especially recommended for use in everyday home use;

It acts on: *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Salmonellae Typhimurium*, *Candida albicans*, *Escherichia coli*, *Bacillus cereus*, *Enterococcus hirae*, *Proteus hauseri*.

METHOD OF USE

DEZISEPT Gel is used by pouring directly from the packaging.

On previously washed and dry hands, apply 3-5 ml of gel and rub it in for 30 seconds. Wait for it to dry. Repeat treatment as needed.

NOTE

Precautionary statements: Keep out of reach of children. Keep away from heat, hot surfaces, sparks, open flames and other sources of ignition. No smoking. Keep container tightly closed. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Store in a well-ventilated place. Keep cool. Disposal of contents / container in accordance with local regulations.

First aid: If swallowed: Rinse mouth with water and drink water (only if the person is conscious). Do not put anything in the mouth of an unconscious person. Seek professional medical attention immediately. Show the doctor the safety data sheet or label.

Dispose of residues of the biocidal product in accordance with the Law on Waste Management (Official Gazette 36/09, 88/10, 14/16 and 95/18-other law), and with used packaging in accordance with the Law on Packaging and Packaging Waste (Official Gazette 36/09 and 95/18-other law).

STORAGE

Protect the product from open flames, heat and direct sunlight. Keep away from food, drink and animal feeding stuffs. Keep away from sources of ignition - no smoking.

Store in a dry, cool and well-ventilated area, in the original packaging.

SHELF LIFE

2 years.

PACKAGING

50ml, 100ml, 300ml

Decision number: 532-01-1322/2020-03 od 21.04.2020.

Manufacturer: PROFARM DOO, Nemanjina 36, 11000 Beograd, tel: 011 4080 108