

GINGIPRO®

Gel with Natural Active Ingredients for Gum Care and Massage

COMPOSITION

Gel contains:

- Essential sage oil 0,6%
- Tormentililae tincture 3%
- Peppermint essential oil 0,6%
- Chamomilae extract 9%

CHARACTERISTICS

The fact of the GingiPro gel is based on the effect of herbal ingredients that are part of its composition:

- * maintains gum health
- * nourishes and favorably affects the gingival mucosa
- * soothes irritated mucous membranes
- * acts astringent
- * acts mildly cool
- * helps prevent gum infections

GingiPro is recommended as an addition to therapy in patients exhibiting obvious symptoms of gums inflammation caused by gingivitis, periodontal disease, accumulated tartar and other.

In dental practice, it is suitable for use in the context of interventions, such as:

- During and period after the placement of mobile prosthetic restorations (total and partial dentures);
- Placement of fixed prosthetic restorations (bridges and crowns);
- After teeth trimming for fixed prosthetic restorations.

METHOD OF USE

Before each application of GingiPro gel, teeth should be washed well with brush and toothpaste. Dispense about 0.5 cm of gel on the finger-tip and apply gentle circular movements to the gel on and around the gums, with gentle circular movements, both on the outside and on the inside, lightly massaging them.

After the gel has been applied, do not take food or liquids for one hour to allow the gel to take action.

GingiPro is used two to three times a day for a minimum of 7 days.

NOTE

Close tube immediately after use.

STORAGE

Do not store above 25°C. Store in original packaging.

SHELF LIFE

2 years.

PACKAGING

A tube with 15g of gel.



11000 Belgrade, Nemanjina 36
Serbia

P165