

FLUOROGAL®

Low Fluoride Concentration Solution for Tooth Protection in Children and Adults

COMPOSITION

1ml of solution contains: sodium fluoride 0,5mg (223 ppmF).

CHARACTERISTICS

FLUOROGAL Solution Red is with strawberry taste. Its daily application enables constant maintenance of fluoride levels on the tooth surfaces, which contributes to the preventive protection of teeth from caries.

Fluorogal Solution is used in adults and children older than 5 years.

METHOD OF USE

Having brush your teeh, pour the Fluorogal solution in amount of one clouser and swirl it around in the mouth for 1 minute. Do not take food or fluid half an hour afterwards.

NOTE

Do not swallow the solution!

Close the bottle immediately after use.

STORAGE

Store at the temperature from 15°C - 30°C.

S2: Keep out of the reach of children.

SHELF LIFE

2 years.

PACKAGING

250ml of FLUOROGAL Solution Red, Strawberry Flavour



d.o.o.

11000 Beograd, Višnjički venac 73D

P019